


Safe n' Dry[®]

Discreet. Dignified. Dependable.

The completely external Safe n' Dry seals gently and securely to the tip of the penis with BioDerm's advanced hydrocolloid adhesive faceplate. Safe n' Dry keeps the skin dry and comfortable—all day, every day.



Safe n' Dry:

- Is completely external
- Keeps users dry and comfortable 24/7
- Fits all male anatomy: large, small, circumcised, uncircumcised and retracted
- Hypoallergenic and latex free
- Seals securely for 24+ hours
- Allows men to intermittent cath with the device in place
- Connects to any standard leg or bedside bag

Part #	Size	Item Per Box	Boxes Per Case
21026	One size	35	12