


External Catheter for Men

A Revolutionary Solution for Male Urinary Management

- Keeps patients dry and comfortable
- One size fits most
- Completely external design
- An occlusive seal around the urinary meatus
- Replaces diapers, condom catheters and medically unnecessary foley catheters
- Suitable for retracted, circumcised and uncircumcised anatomies
- Avoids moisture related skin breakdown




Men's Liberty™
with drip chamber



FreeDerm™ Adhesive Remover
Instantly removes hydrocolloid without irritation.



BioPlus™
Barrier Prep Wipe



CathGrip® UniGrip™
The recommended device to anchor Men's Liberty Acute

Decrease Infection Rates Today!

- Indwelling Catheters increase risk of infection (80% of UTI's are associated with use of an indwelling catheter which lengthen hospital stays and increase costs)
- Condom Catheters can result in leakage and skin breakdown
- **Men's Liberty Delivers Superior Performance with leak protection, Ease of Use and Patient Comfort**

Wear Time: The first application for each patient will last shorter than the subsequent applications of the product. You can expect the first application to last 6-12 hours and subsequent applications to last 24-48 hrs +. This is due to the hydrocolloid material that absorbs underlying moisture in the skin.

For more information please contact your sales representative.

| Item Description | Medline Item # | Unit of Measure | Per Box |
|--------------------------------|----------------|-----------------|---------|
| Men's Liberty Acute | BIO210268 | 1 Box | 35 |
| Men's Liberty Acute w/CathGrip | BIO2102610 | 1 Box | 35 |
| Men's Liberty - w/drip chamber | BOM23046BX | 1 Box | 35 |
| CathGrip Large Double Strap | BOM51300NSBX | 1 Box | 30 |
| CathGrip Medium Double Strap | BOM51301NSBX | 1 Box | 30 |
| CathGrip Small Single Strap | BOM51250NSBX | 1 Box | 30 |
| FreeDerm (1mL ampoules) | BOM52203BX | 1 Box | 30 |
| BioPlus+ Wipes | BOM52300BX | 1 Box | 25 |

Removal: Remove and replace MLA when the adhesive turns a white color, when it is loose, or when it has been soiled. We recommend using Freederm adhesive remover to remove Men's Liberty Acute.



Manufactured by BioDerm, Inc.
12320 73rd Court North, Largo, FL 33773
Email: CustomerCare@BioDermInc.com
Web: BioDermInc.com

© 2018 BioDerm, Inc. All rights reserved. BioDerm® and CathGrip® are registered trademarks of BioDerm, Inc. BioPlus+™, FreeDerm™, Men's Liberty™ and Men's Liberty Acute™ are trademarks of BioDerm, Inc. All other trademarks acknowledged. U.S. patents #8,551,062, #D767,761. CathGrip® and Men's Liberty Acute™ U.S. patent pending.



Medline is the largest privately held manufacturer and distributor of health care supplies in the United States.

1-800-MEDLINE (633-5463) | www.medline.com