



**PRODUCT CATALOG**




## Innovative solutions in urology, skin protection and securement

BioDerm® is one of the nation's leading manufacturers of proprietary hydrocolloid products for securement, urinary management and skin protection. The company challenges accepted inferior standards of care by creating products that reduce infection rates, add comfort and reliability, and vastly improve quality of life.

### UROLOGY

Men's Liberty™ and Men's Liberty Acute™ are revolutionary alternatives to diapers, condom catheters and when medically unnecessary, Foley catheters.




### SKIN PROTECTION

Ensure proper product application and removal with the use of BioPlus+™ Barrier Film Wipes and FreeDerm® Adhesive Remover.



### SECUREMENT

The CathGrip® family of universal medical tubing securement options provides clinicians with the ability to manage a wide variety of catheters and tubes ranging in size from 6-42 Fr.





## Say Goodbye To Adult Diapers Once And For All

Staying dry and odor-free shouldn't even be a question. The Men's Liberty line of products are non-invasive, external urine collection devices that can be used at home and on-the-go, so you can forget about bladder leakage and get back to living life on your own terms.


- Replaces diapers, condom catheters and when medically unnecessary, Foley catheters
- Has a completely external design that fits most male anatomy
- Keeps men dry and comfortable for up to 24 hours
- 8 oz collection pouch can be emptied easily throughout the day
- Connects directly to a leg or bed bag
- Covered by most insurance policies (Standard deductibles apply)




### How Does It Work?


Men's Liberty & Men's Liberty Acute are made from hydrocolloid, an adhesive that creates a secure, skin friendly seal around the urinary opening. Men's Liberty is designed for use in the outpatient setting and includes an integral collection chamber that holds up to 8 fl. oz (250 mL) of urine. Patients can fill the collection chamber, which can then be discreetly emptied into a toilet through the end plug located on the bottom of the product, whenever they get the opportunity.

#### Position & Apply



#### Secure & Seal






SKU# 23046

## Men's Liberty™

### For Outpatient Use

Men's Liberty is the completely external catheter that features a micro-vent and is applied just once a day. This discreet, one-size-fits-most option manages men's urinary incontinence in a safe, healthy way. More than 5 million units have been sold with zero attributable adverse events confirmed, including UTIs. Men's Liberty is covered by Medicare and over 3,000 insurance plans in the US.



SKU# 23046C

## Men's Liberty™ Classic

### For Outpatient Use

Our Men's Liberty Classic is the same product as Men's Liberty with a larger vent valve. The larger valve helps patients who have difficulty placing the vent up during application. innovative alternative to traditional male external catheters. It is also a completely external, discreet, one-size-fits-most option that manages men's urinary incontinence in a safe, sanitary way.



SKU# 21026

## Men's Liberty™ Acute

### For Acute Care Use


Men's Liberty Acute is a male external catheter designed for institutional and hospital use. It is completely external and seals gently and securely to the tip of the penis with BioDerm's proprietary hydrocolloid adhesive. Men's Liberty Acute is exclusively distributed in the US by Medline Industries. Contact Medline at [customerservice@medline.com](mailto:customerservice@medline.com) or 1-800-Medline.



## Penile Clamp for Men

The KindKlamp, penile clamp, is designed to temporarily stop the flow of urine through the urethra when applying Men's Liberty. This simple accessory helps men who are dribbling or leaking and are having trouble applying the Men's Liberty device.

- Easy use clasp – can be closed and released with one hand!
- Curved upper arm for improved comfort
- Long lasting hinge designed for multiple cycles
- Designed to fit most male anatomy
- The plastic design allows for use during x-rays and diagnostic imaging



SKU# 20003

### KindKlamp™

#### For Outpatient Use

The KindKlamp is made from durable plastic and ultra soft foam for added comfort. It is made without wires or metal. This ensures that the KindKlamp is skin friendly and can be wiped with a damp cloth without the potential of rust or corrosion. One clamp is covered by insurance every 90 days with a doctor's order\*.

\*Standard deductibles and co-pays apply



**1000 mL Leg Bag**  
**SKU# 90000**

## Leg Bag

The Leg Bag can be used for extra capacity while using the Men's Liberty. The universal tube directly connects to most catheter systems. To use, simply unplug the bottom of the Men's Liberty pouch and connect to the Leg Bag.

- Kink-resistant tubing
- All fluid pathways are sterile
- Convenient drain spout
- Straps included
- Anti-reflux valve




**2000 mL Bed Bag**  
**SKU# 90001**

## Bed Bag

The Bed Bag can be used for extra capacity while using the Men's Liberty during the night. The universal tube directly connects to most catheter systems. To use, simply unplug the bottom of the Men's Liberty pouch and connect to the Bed Bag. Two bags are covered by insurance every month with a doctor's order\*.

- Kink-resistant tubing
- All fluid pathways are sterile
- Convenient drain spout
- Hook included to connect to bed railing
- Anti-reflux valve

# SKIN PROTECTION



## Faceplate Strip

Faceplate Strip prevents skin injuries and irritation due to face masks, glasses, and face shields. Face masks, glasses and face shields are a necessity and are being worn more frequently and for longer periods of time. No matter how well your face PPE fits, constant wear time can cause skin injuries and irritation. The wing design will cover the bridge of your nose, cheek bones and forehead for up to 24 hours.

- Made in the USA and are skin friendly
- Can be worn for up to 24 hours



## Bio Plus+<sup>TM</sup> Barrier Film Wipe

BioPlus+ improves hydrocolloid adhesion by up to 50%. A BioPlus+ wipe is included in most BioDerm products and is also available as a stand alone product.

- No-sting
- Forms transparent, protective barrier film
- Dries faster than major competitors in lab tests
- Increases skin adhesion more than major competitors in lab tests





## Adhesive Remover

Specially formulated to instantly remove hydrocolloid, tape and bandages.

Instantly removes hydrocolloid without irritation and evaporates rapidly with no residue. Allows immediate re-application of an adhesive device without washing. Reduces skin trauma from tape and dressing removal including infection risk, inflammatory damage and edema. Available as 1 oz or 3 oz pump spray, single-use wipes or ampoules.

- Safe and gentle for all patients
- No-sting
- Instantly removes hydrocolloid
- No harsh chemicals
- Oil and lotion free

**Dries in 15 seconds**  
allowing for immediate reapplication



# SECUREMENT

---



## The CathGrip Family of Universal Securement Options

The CathGrip line provides clinicians with the ability to manage a wide variety of catheters and tubes ranging in size from 6-42 Fr. Each product is based on BioDerm's proprietary hydrocolloid with features that solve a wide variety of challenges. When properly applied, these devices provide a comfortable fit and secure hold for up to seven days. The soft straps are durable for repeat positioning, perfect for continual ease of use. CathGrip can meet all your securement needs, making it the new standard of care sought after by physicians, nurses and purchasing teams.

**The CathGrip family of products is exclusively distributed in the US by Medline Industries. Contact Medline at [customerservice@medline.com](mailto:customerservice@medline.com) or 1-800-MEDLINE.**




**Not made with  
natural latex**



**Up to 7-day  
wear time**




**Fits 6-42 Fr.  
medical tubing**





# Single Strap & Double Strap

CathGrip is a truly universal system that securely holds tubes from 6-42 Fr. Available in single or double straps, CathGrip provides all over securement with a skin-friendly hydrocolloid adhesive to prevent skin injuries. Soft straps are durable for repeat positioning.




- Available in single or dual straps in the medium and large sizes
- Hydrocolloid adhesive oval
- Not made with natural rubber latex
- Soft straps with no-slip grip
- Simple repeat positioning
- Fits all medical tubing from 6-42 Fr
- Durable repeat positioning
- Extended wear time, 2-3 times longer than competitors

## Sterile & non-sterile available



## Low Profile Single Strap & Double Strap



A new addition to the CathGrip family. The Low Profile CathGrip has all of the qualities of the existing Single and Double Strap but with the addition of low profile straps. The low profile straps hug the body reducing the risk of snagging for a more comfortable wear. Sterile and non-sterile available.



## Oval

Designed to prevent pressure injuries and skin tearing with low-profile proprietary hydrocolloid securement, especially important when dealing with tender skin. The ovals are sold in pairs. The low profile design allows for application in difficult areas. Sterile and non-sterile available.

Uses listed above are not all inclusive. May vary based on tube location, infection prevention protocols and presence of exuding wounds.




# Single Strap with Wing Seal

## Double Strap with Wing Seal

Extends the security of our hydrocolloid base to support insertion site securement.

Hydrocolloid wings form an occlusive seal that securely holds any tube from 6-42 Fr. and can help reduce dislodgement and minimize the risk of infection. Only available as sterile.




Uses listed above are not all inclusive. May vary based on tube location, infection prevention protocols and presence of exuding wounds.

- Available in single or double straps
- Hydrocolloid adhesive oval
- Not made with natural rubber latex
- Soft straps with no-slip grip
- Simple repeat positioning
- Fits all catheters from 6-42 Fr
- Durable repeat positioning
- Extended wear time, 2-3 times longer than competitors



**Single Strap with Wing Seal**  
SKU# 51100





**Double Strap with Wing Seal**  
SKU# 51303

# PEG Protect™

Designed specifically for use with PEG tubes, PEG Protect brings our skin friendly hydrocolloid securement to long and short term PEG patients. Only available as sterile.

- Allows for regular tube rotation for traditional and low profile PEG tubes
- Reduces the risk of dislodgement and tube migration
- Provides a skin friendly hydrocolloid barrier against bodily fluid, bruising and excoriation
- Reduces potential for tube kinking
- Two piece design allows for variable securement location related to the PEG tube insertion site

**Secures both standard and low-profile PEG tubes:  
allowing regular tube rotation for traditional and  
low-profile PEG tubes, Nephrostomy Chest Tube,  
PEG Tube and J Tube**



SKU# 51501



**The CathGrip family of products is exclusively distributed in the US by Medline Industries. Contact Medline at [customerservice@medline.com](mailto:customerservice@medline.com) or 1-800-MEDLINE.**

# Now That You Know Our Products

BioDerm provides innovative solutions for customers around the world. Safer, easy to use and cost effective, our proprietary hydrocolloid, skin prep, adhesive remover and soft strap technology is the next evolution in patient care. Join the future of healthcare technology today!

## OUTPATIENT

BioDerm's products are sold through our customer service center, based in Largo, Florida. Our friendly customer care team is trained to understand the conditions associated with urinary incontinence and are waiting for your call.

## ACUTE CARE

BioDerm's CathGrip and Men's Liberty Acute products are sold exclusively through Medline who targets acute care, long-term care and other surgical centers in the U.S. The CathGrip line of products is transforming care in critical care units, interventional radiology, medical/surgical units, wound care clinics and more with solutions that help reduce the incidence of medical device related pressure injuries and catheters becoming dislodged.

## HEALTHCARE PROVIDERS

BioDerm works seamlessly with many HCP offices to help patients obtain insurance coverage for Men's Liberty. Contact us today for a HCP Provider Kit.

## INTERNATIONAL

BioDerm's products are currently available for sale in over 10 countries through a network of independent distribution partners in each country.

## PARTNERSHIPS

BioDerm products are available for private label collaboration for target markets and partnerships. Contact our care team today.





## Contact Us

### **For Our U.S. Customers**

Call 1-800-814-3174

Email [CustomerCare@BioDermInc.com](mailto:CustomerCare@BioDermInc.com)

Fax completed LMNs to (888) 920-9370

### **For International Inquiries**

Please contact [CustomerCare@BioDermInc.com](mailto:CustomerCare@BioDermInc.com)

# ADDITIONAL INFORMATION

## Regulatory Information

BioDerm is registered with the US Food and Drug Administration and other international regulatory bodies. Copies of registration documents are available upon request or online at <http://bioderminc.com/educational-resources/product-documentation/>

| Product Category   | Role                                |
|--|-------------------------------------|
| Tape and bandage, adhesive – BioPlus+ Barrier Film Wipe  | Specification Developer             |
| Device, intravascular catheter securement – CathGrip Double Strap; Cathgrip Double Strap with Wing Seal; CathGrip Oval; CathGrip PEG Protect; CathGrip Single Strap; CathGrip Single Strap with Wing Seal; Tube Securement Device; Low Profile CathGrip Single Strap; Low Profile CathGrip Double Strap; Medium LVAD Securement Device; Large LVAD Securement Device | Contract Manufacturer; Manufacturer |
| Solvent, adhesive tape – FreeDerm Adhesive Remover; FreeDerm Bottles; FreeDerm Ampoules; FreeDerm Wipes  | Specification Developer             |
| Device, paste-on for incontinence, non-sterile - Faceplate Seals; Men's Liberty Acute with CathGrip; Men's Liberty Acute with Cathgrip and Towel; Men's Liberty Acute without CathGrip; Men's Liberty; Men's Liberty Acute; Men's Liberty Classic  | Contract Manufacturer; Manufacturer |
| Clamp, penile - BioDerm Penis Clamp; KindKlump   | Specification Developer             |
| Bandage, elastic - Faceplate Strip   | Manufacturer                        |



\*does not apply to leg bag, bed bag and sterile products.

### US Food and Drug Administration

Establishment Registration Number: 1063299  
 BioDerm, Inc., 12320 73rd Court North,  
 Largo, FL 33773  
 Status: Active

### Authorized Representative: Europe

Advena Ltd  
 Tower Business Centre, 2nd Flr,  
 Tower Street, Swatar, BKR 4013 Malta  
[www.advenamedical.com](http://www.advenamedical.com)

### British Standards Institute

Certificate of Registration:  
 Quality Management  
 System ISO-13485:2016  
 Certificate #: FM 565097  
 Originally Registered: 07/06/2011

| Product Name                                | HCPSC Code                             | Catalog/SKU Number* | Size          | Units Per Box | Boxes Per Case | Expiration Date<br>(from the date of manufacture)* |
|---|--|---------------------|---------------|---------------|----------------|--|
| <b>UROLOGY</b>                              |  |                     |               |               |                |  |
| Men's Liberty                               | A4326                                  | 23046               | One Size      | 35            | 12             | 30 months*   |
| Men's Liberty Classic                       | A4326                                  | 23046C              | One Size      | 35            | 12             | 30 months*   |
| Men's Liberty Acute                         | NA                                     | 21026               | One Size      | 30            | 12             | 30 months*   |
| KindKlamp                                   | A4356                                  | 20003               | One Size      | 12            | n/a            | n/a  |
| Drainage Bag – Leg                          | A4358                                  | 90000               | 1000 mL       | 50            | 4              | 60 months  |
| Drainage Bag – Bedside                      | A4357                                  | 90001               | 2000 mL       | 50            | 4              | 60 months  |
| Faceplate Seals (Pack of 10)                | NA                                     | 20002               | One Size      | 1 Pack        | n/a            | 48 months  |
| <b>SKIN PROTECTION</b>                      |  |                     |               |               |                |  |
| BioPlus+ Wipe                               | A5120                                  | 52300               | One Size      | 25            | 12             | 30 months  |
| FreeDerm Spray (1oz)                        | A4455                                  | 52200               | 1 oz.         | 36            | 4              | 36 months  |
| FreeDerm Spray (3oz)                        | A4455                                  | 52202               | 3 oz.         | 63            | 1              | 36 months  |
| FreeDerm Ampoule                            | A4456                                  | 52203               | 1 mL          | 30            | 36             | 36 months  |
| FreeDerm Wipe (6mL)                         | A4455                                  | 90037               | 6 mL          | 30            | 1              | 36 months  |
| Faceplate Strip                             | NA                                     | 20002NP             | 4.5" x 1"     | 10            | 100            | 48 months  |
| <b>SECUREMENT - CATHGRIP</b>                |  |                     |               |               |                |  |
| Double Strap Large (2 Straps)               | A4333 for Urinary                      | 51300               | 4" x 4.75"    | 30            | 3              | 30 months*   |
| Double Strap Large (2 Straps, Non-Sterile)  |  | 51300NS             | 4" x 4.75"    | 30            | 3              | 30 months*   |
| Single Strap Large (1 Strap)                |  | 51350               | 4" x 4.75"    | 30            | 3              | 30 months*   |
| Single Strap Large (1 Strap, Non-Sterile)   |  | 51350NS             | 4" x 4.75"    | 30            | 3              | 30 months*   |
| Double Strap Medium (2 Straps)              |  | 51301               | 3" x 3.5"     | 30            | 3              | 30 months*   |
| Double Strap Medium (2 Straps, Non-Sterile) |  | 51301NS             | 3" x 3.5"     | 30            | 3              | 30 months*   |
| Single Strap Small (1 Strap)                |  | 51250               | 3" x 3.5"     | 30            | 3              | 30 months*   |
| Single Strap Small (1 Strap, Non-Sterile)   |  | 51250NS             | 3" x 3.5"     | 30            | 3              | 30 months*   |
| Low Profile Single Strap (1 Strap)          |  | 51250LP             | 3" x 3.5"     | 30            | 3              | 30 months*   |
| Low Profile Double Strap (2 Straps)         |  | 51300LP             | 4" x 4.75"    | 30            | 3              | 30 months*   |
| Double Strap with Wing Seal                 | A5200 for Percutaneous                 | 51303               | 2.94" x 4.75" | 30            | 3              | 30 months*   |
| Single Strap with Wing Seal                 |  | 51100               | 2.2" x 3.56"  | 30            | 3              | 30 months*   |
| PEG Protect                                 |  | 51501               | 3" x 3.5"     | 30            | 3              | 30 months*   |
| Oval Large                                  | Varies Based on Type of Tubing Secured | 54003               | 3" x 3.5"     | 30            | 3              | 30 months*   |
| Oval Large (Non-Sterile)                    |  | 54003NS             | 3" x 3.5"     | 30            | 3              | 30 months*   |
| Oval Medium                                 |  | 53002               | 1.48" x 2.4"  | 30            | 3              | 30 months*   |
| Oval Medium (Non-Sterile)                   |  | 53002NS             | 1.48" x 2.4"  | 30            | 3              | 30 months*   |
| Oval Small                                  |  | 54002               | .93" x 1.5"   | 30            | 3              | 30 months*   |
| Oval Small (Non-Sterile)                    |  | 54002NS             | .93" x 1.5"   | 30            | 3              | 30 months*   |

\*Expiration date is based on the earliest expiration date of the components. Expiration date indicates the expiration of the included BioPlus+ wipe. Products are not shipped within 1 year of expiration.

CathGrip, Men's Liberty and Men's Liberty Acute ship without alcohol wipes in Canada. Product availability is dependent on the regulations of your country. Not all products are available in all countries.



12320 73rd Court North, Largo, FL 33773

US Sales: 800-814-3174

International Sales: +01 727-507-7655

Support: 800-864-6755

Email: [CustomerCare@BioDermInc.com](mailto:CustomerCare@BioDermInc.com)

Web: [BioDermInc.com](http://BioDermInc.com)

© 2021 BioDerm, Inc. All rights reserved. BioDerm®, FreeDerm® and CathGrip® are registered trademarks of BioDerm, Inc. BioPlus+™, KindKlamp™, PEG Protect™, Men's Liberty™ and Men's Liberty™ Acute are trademarks of BioDerm, Inc. Other trademarks acknowledged. U.S. patents #D7,896,857, #8,551,062, #10,022,261, #10,485,955, #767,761, #D890,338, #D890,337, #10,369,332. International Patents: (Canada) #3,058,756, #30,595,573 (UK and NE) 2501429, (UK) 1827530 and (Japan) # 4939422. Additional patents pending.

LT 93943 F